

B-FORCE 4.0

FOR A DOUBLE
SIMULTANEOUS ACTION
ON MUSCLES AND FAT



BALDAN
G R O U P

GET IN SHAPE
IN NO TIME

NOT ONLY BEAUTY, THE IMPORTANCE OF A HEALTHY BODY

To maintain a healthy lifestyle and enjoy immediate and long-term benefits, the muscular system must be constantly trained.

TONING THE BODY OFFERS MULTIPLE BENEFITS FROM AN AESTHETIC AND, ABOVE ALL, PHYSICAL POINT OF VIEW.

Muscle strength training brings a general sense of well-being, significantly improves quality of life, helps control anxiety and stress, reduces the risk of chronic diseases, and stabilizes bones and joints.

Finding the right motivation is a problem for many people who, most of the time, prefer to spend their free time in other ways to escape the daily stress, thus putting aside their physical well-being.

Baldan Group once again breaks the barriers of wellness and erases the negative aspects of exercise with B-FORCE, the technology designed specifically for those who want to stay in shape or regain their ideal figure in just a few weeks

Working on your body is important

A new training concept for the **WELL-BEING OF THE BODY** at 360 degrees

Research has led to the creation of B-FORCE, the first electrostimulator which, thanks to a unique combination of various waveforms, enables 5 different types of contractions to increase muscle tone and reduce localized fat deposits.

B-FORCE is the first technology that allows variable contractions, unlike other similar devices that offer only one type of contraction, more or less intense.

The contractions triggered by electrostimulation stimulates muscular energy which, as it continues to be stimulated, increases the volume of the treated muscle, attacks fat and increases metabolism.





Ideal for anyone wanting to work on localized fat.

ELECTROSTIMULATIONS

Electrostimulation sends direct signals to the muscle fibers, which are strengthened: the complete contractions involve the muscle in depth, resulting in a simultaneous action on all the muscle tissue in the treated area.

With B-FORCE, excess fat is reduced and, at the same time, the treated area is firmed and sculpted.

THIS LATEST-GENERATION TECHNOLOGY WORKS SAFELY, QUICKLY AND PAINLESSLY ON REMODELLING AND TONING OF THE BODY.

THE DESIGN

Completely MADE IN ITALY, the designer of the project is Simone Micheli, one of the most famous architects in the world, called to express concretely the innovative character of the device: to communicate solidity, Micheli used the Baidur 100 material, patented by Bayer.



Thanks to its dual action, B-FORCE acts on muscles and at the same time combats fat.

A B - FORCE session lasts 20 to 30 minutes!

B-FORCE 4.0

THE SOFTWARE

Thanks to the adoption of a state-of-the-art computer intelligence system, the practitioner can constantly monitor the progress of treated clients and missed sessions, so as to stay up-to-date.

Personal evaluation forms can be created on the basis of automatic indications obtained from an ideal protocol, supplied directly by Baldan Group, to be adopted for each type of client according to pre-set forms based on gender, area to be treated and type of imperfection.



**MONITORING OF
CUSTOMER PROGRESS**



**MADE-TO-MEASURE
DATA SHEET**

B-FORCE IS A SMART TECHNOLOGY, THANKS TO THE INTERNET CONNECTION INTEGRATED IN THE DEVICE.

THIS ALLOWS CONSTANT DIALOGUE WITH BALDAN GROUP'S EXPERTS, WHO CAN AUTOMATICALLY REPORT MALFUNCTIONS OR USAGE PROBLEMS, INSTANTLY ADD NEW PROGRAMMES AND CARRY OUT REMOTE ASSISTANCE


An innovative aesthetic technology that works on reshaping the body.

The SIMPLIFICATION OF SOFTWARE for enhanced treatment security

The simplification of the software ensures greater treatment safety, thanks also to the development of automatic protocols offering pre-set parameters and working times divided into zones, to ensure more efficient practices and targeted treatment of the most likely imperfections.

Once the data has been entered, all useful information for the operator such as the handpiece to be used, the pulse frequencies and the treatment time will appear on the 15-inch touch screen.

IT IS ESSENTIAL THAT DATA IS UPDATED SYSTEMATICALLY, SESSION AFTER SESSION.



Acts on Abdomen, Thighs,
Buttocks, Arms, Calves and Hips

HANDLES

The technology acts on different areas of the body through its 4 handpieces:

2 HANDLES for treating large and medium-sized areas such as the ABDOMEN, THIGHS, GLUTES

2 anatomically shaped HANDLES to work on smaller, curved areas, such as ARMS, CALVES AND HIPS.

B-FORCE 4.0





The programme reaches deep into the muscles, even those that are usually difficult to strengthen during traditional training.

An air-cooling system is built into each applicator so that the handpieces are lighter for both the operator and the client undergoing treatment.

B-FORCE 4.0

No Pain, no Damage to the surface,
no Recovery time.




A B-FORCE SESSION is
FAST AND HARMLESS
for the customer, at HIGH
PROFITABILITY for the
entrepreneur

It is a highly profitable activity for the practitioner, as there is no consumption of materials and products and, moreover, it does not require the continuous presence of the practitioner during the course of the session.


THE MULTIFUNCTIONAL TECHNOLOGY THAT RESHAPES THE BODY IN A PERSONALISED MANNER

B-FORCE is capable of solving various problems, just like an equipment room in a gym where the customer can choose which workout to do according to their needs and requirements.


The machine incorporates innovative technical aspects to ensure certain, reliable and effective results, guaranteeing greater efficiency and optimum performance.



**REDUCTION OF
ADIPOSE TISSUE**



**MEASURABLE TONING OF
TREATABLE MUSCLE**



**REDUCTION IN
CIRCUMFERENCE**

THE 7 ADVANTAGES OF B - FORCE

/01

DRAINAGE and MICROCIRCULATION

The contractions are stepped and this causes a vibratory contraction that acts on the lymphatic vessels, promoting drainage and microcirculation.

This programme is suitable for treating adiposity and localized cellulite.

/02

SCULPTING

This programme involves deep contractions and de-contractions that force the muscle to perform intense work. The programme is suitable for sculpting the body and reducing fat at the same time.

/03

MUSCLE TONING

In this programme, contraction takes place more slowly, the muscle remains contracted for a longer period of time and is quickly released. This type of activity increases muscle tone and prepares the ground for future muscle development. This programme is recommended in the case of poor muscle tone.

/04

VOLUME MUSCULAR

Contraction occurs progressively, the muscle remains contracted for a medium length of time and de-contraction occurs progressively. This type of activity produces an increase in white muscle fibrils and, consequently, an increase in muscle volume. It is suitable for people who have few muscles and want to develop their mass.

/05

STRENGTH

Contraction and de-contraction follow each other with little pause and without any kind of jolt. The muscle is practically subjected to high levels of intensity. Suitable for athletic individuals who have good muscle mass and want to increase strength and endurance.

/06

CIRCUIT 1

Work protocol that, session after session, alternates all the programmes listed above. Suitable for subjects who have several problems to solve at the same time, circuit indicated mainly for slimming.

/07

CIRCUIT 2

For those who want a complete 'equipment room' within the institute to improve their muscles as a first objective. Suitable for athletes or those who love physical activity.

PARAMETERS

It is important to act directly on the specific needs of the customer to guarantee the best results, which is why all programmes can be customised according to the following parameters:

- ✓ WAVE SHAPE
- ✓ PULSE WIDTH
- ✓ INTENSITY
- ✓ FREQUENCY
- ✓ RELAXATION TIME
- ✓ CONTRACTIONS TIME

TREATMENT RECOMMENDED FOR:

- ☐ INCREASING MUSCLE MASS AND TONE
- ☐ REDUCING ADIPOSE TISSUE
- ☐ TONIFYING
- ☐ COMBATING THE MOST STUBBORN IMPERFECTIONS
- ☐ STRENGTHENING BONES AND JOINTS
- ☐ IMPROVING POSTURE

PROTOCOLS OF USE

TIMING: the average time required for a treatment with B-FORCE is about 20/30 MINUTES. per area.

EXPERIENCE: most clients relax comfortably during the treatment.

REGIME: Treatments are generally scheduled on the basis of 2 sessions per week. Most clients will sustain 8 to 10 treatments to achieve the desired result.



B-FORCE 4.0

BALDAN
G R O U P
www.baldangroup.com